



GOVERNMENT OF WEST BENGAL
OFFICE OF THE DIRECTOR
INSTITUTE OF POST GRADUATE MEDICAL EDUCATION & RESEARCH
244, A.J.C BOSE ROAD, KOLKATA- 700020

Memo No. Inst/ 8639

Dated: 07.06.18

NOTICE

Sealed quotations are invited from bonafide/interested firms with superscription "Quotation for project ICMR, Govt. of India entitled 'Identification of clinico-immunological determinants.....Leishmaniasis (PKDL)'. Ref: 6/9-7(151)/2017-ECDII. The quotation should reach the office of the undersigned by 18/6/18.....at 2 pm. The quotation should accompany with current valid trade License, PAN No., GST Registration Certificate proof.

Acceptance of the lowest rate is not obligatory for the undersigned that reserves the right to accept or reject any offer without assigning any reason whatsoever. Rate must be submitted with company name & cat no. after calculating GST or other Taxes, if any according to the following format.

As per letter of Principal Investigator of the project these items need to be purchased as per specifications mentioned, the quotation is placed herewith. The quotation may be opened on or after 19/6/18..... onward, bonafide/interested firm's representative may contact Store In-charge (Projects).

1. FACS Lysing solution: 100 ml

- Should be intended for direct immunofluorescence staining of human peripheral blood cells
- Should be compatible for both lyse/wash and lyse/no-wash procedures

2. RBC lysis solution: 100 ml

- Should be ammonium chloride-based lysing reagent
- Should not contain a fixative agent, so leukocytes remain viable after red blood cell lysis

3. FACS Flow Sheath Fluid: 20L

- Should contain sodium fluoride (NaF) and be compatible with BD FACS Verse Flow cytometer

4. FACSuite CS&T Beads 50Tests (RUO):

- CS&T research beads should consist equal quantities of 3- μ m bright, 3- μ m mid, and 2- μ m dim polystyrene beads in phosphate buffered saline (PBS) with bovine serum albumin (BSA) and 0.1% sodium azide
- Should be compatible with BD FACS Verse Flow cytometer

5. Fixation/Permeabilization Solution: 250 Tests

- Should have a property for simultaneous fixation and permeabilization of cells prior to intracellular cytokine staining and be compatible with BD FACS Verse Flow cytometer

6. Perm/Wash Buffer:

- Should contain Fetal Bovine Serum (FBS) and saponin and be compatible with BD FACS Verse Flow cytometer

7. Specification for Microcentrifuge

- Max speed 13,000 rpm and time control (pulse <30 seconds)
- Rotors should have capacity of 6 (x 1.5 ml tubes) and 12 (X 0.5 ml tubes)


Director,
IPGMER, Kolkata SR.
Director
I. P. G. M. E. & R.
Kolkata-700020